

## Overview

Alleviate adulteration and tampering concerns by upgrading to the iCup device. The iCup® Drug Screening Device simplifies the urine collection process while providing fast and reliable test results. Numerous devices are available, covering up to 10 common drugs of abuse. Adulteration test strips also available. Provides instant test results for a wide range of illicit drugs. The self-contained cup is ideal for sending presumptive positive specimens to the lab for confirmation.

- Test Results Can Be Photocopied
- Eliminates the Need for Restrooms
- Difficult to Adulterate
- Easily Transport to Labs
- Drug Test results are Displayed in Minutes!
- Do Not Place Anything in the Mouth Including Food Drink or Gum prior to the Test
- Negative and Positive Results are Clear and Easy to Read

<b>TEST</b>	<b>CODE</b>
Amphetamine	AMP
Barbiturates	BAR
Benzodiazepines	BZO
Buprenorphine	BUP
Cocaine	COC
Synthetic Cannabinoids	K2/SPICE
Marijuana	THC
Methadone	MTD
Methamphetamine	M-AMP
Morphine	MOR
Opiates	OPI
Alcohol	ETG
Phencyclidine	ETG
Tricyclic Antidepressants	TCA
Ecstasy	MDMA
Oxycodone	OXY
Propoxyphene	PPX
Fentanyl	FTY
Tramadol	TRAM

Minimum 12-month shelf life.

## Attributes

Option AMP, BZO, COC, THC, M-AMP, OPI, BAR, BZO, THC, M-AMP, OPI, MDMA, OXY,  
AMP, BAR, BZO, COC, THC, MTD, M-AMP, OPI, ETG, TCA, AMP, BAR, BUP, BZO,  
COC , THC, MTD, M-AMP, OPI, ETG, TCA, OXY,PPX

Packed 25/Cs

Product Unit Type Case

Images

